


LASIK CANNULA

- For irrigation under the flap
- Designed for direct removal of flap debris
- Prevents reflux of unwanted fluid under the flap



REF #	DESCRIPTION
4036L	Lin LASIK Cannula 25g x 7/8", 9mm flattened angle, two 0.3mm ports
4036LB	Buratto LASIK Cannula 25g, 11mm angle, one front and two side ports
4036LR	Lindstrom LASIK Irrigator 25g, 45° angle, 5/8" radius Three 0.3mm ports
4036LS	Slade LASIK Irrigator 26g x 1" 9mm formed extension with spatulated tip.
4123	LASIK Cannula 23g x 7/8", Flattened Bend to Tip
4225	LASIK Cannula 25g x 7/8", Formed 7mm from Open End Six 0.15mm side ports at 2mm, 4mm and 6mm from end